

**Molecular  
Devices**

Together through life sciences.

## MetaXpress® Software: *Angiogenesis Module*



Together through life sciences.

©2014 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.


 **Molecular  
Devices**

# Angiogenesis

- Identifies and measures tubes (long thin objects) and nodes (connecting points between tubes) in a single wavelength
- In addition to tube formation assays, can be used for to measure neurite outgrowth in assays where the cell bodies are indistinct or outside the field of view




# Module settings




- Source image
- Select the wavelength for the **tubes** and **nodes**

# Module settings




- Display result image
- Leave “Display result image” deselected (this is generally only used when journaling)



# Module settings



- Parameters
- Set the **Approximate min width** and **Approximate max width** for the range of tubes that you want to detect
- The width is the distance across a tube (in  $\mu\text{m}$ )
- The region tools can be used to measure widths
- Much smaller tubes will be ignored
- Much larger tubes will be identified as nodes





# Module settings




- Parameters
- The **intensity above local background** is used for finding the tubes and nodes
- This value is a minimum and should be set slightly lower than the difference in intensity between the dimmest part of the tube (typically the thinnest part) and its local background
- Draw a line across a cell into the background and use Measure > Linescan to determine intensity values; or simply mouse over the cell and the background and view the intensity values

# Module settings



- **Segment detection**
- **Watershed lines (original)** uses the legacy segmentation method that was used in MetaMorph prior to version 7.6.4.0 and MetaXpress prior to version 3.1.0.69
- **Skeletonization** uses a new segmentation method, which may provide more accurate segmentation for some data sets that have thicker and sparser tubes

# Regions for Measurements





# Summary Data (site-by-site measurements)

- ✓ Total Tube Length
- ✓ Mean Tube Length
- ✓ Total Tube Area
- ✓ Mean Tube Area
- ✓ Tube % Area Covered
- ✓ Average Tube Thickness
- ✓ Segments
- ✓ Branch Points
- ✓ Nodes
- ✓ Total Node Area
- ✓ Mean Node Area
- ✓ Node % Area Covered
- ✓ Connected Sets
- ✓ Tube Length Per Set

- **Total Tube Length:** Total microns of tube length (excluding nodes)
- **Mean Tube Length:** Total tube length divided by the number of segments
- **Total Tube Area:** Total square microns of tube area (excluding nodes)
- **Mean Tube Area:** Total square microns of tube area divided by the number of segments

# Summary Data (site-by-site measurements)

- ✓ Total Tube Length
- ✓ Mean Tube Length
- ✓ Total Tube Area
- ✓ Mean Tube Area
- ✓ Tube % Area Covered
- ✓ Average Tube Thickness
- ✓ Segments
- ✓ Branch Points
- ✓ Nodes
- ✓ Total Node Area
- ✓ Mean Node Area
- ✓ Node % Area Covered
- ✓ Connected Sets
- ✓ Tube Length Per Set

- **Tube % Area Covered:** Total tube area (excluding nodes) divided by total image area (width times height), times 100
- **Average Tube Thickness:** Average thickness of tubes computed as total area (excluding nodes) divided by total length (excluding nodes), and indicated in microns
- **Segments:** Total number of tube segments connecting branch points and/or ends
- **Branch points:** Total number of junctions connecting segments (nodes are not considered branches)

# Summary Data (site-by-site measurements)

- ✓ Total Tube Length
- ✓ Mean Tube Length
- ✓ Total Tube Area
- ✓ Mean Tube Area
- ✓ Tube % Area Covered
- ✓ Average Tube Thickness
- ✓ Segments
- ✓ Branch Points
- ✓ Nodes
- ✓ Total Node Area
- ✓ Mean Node Area
- ✓ Node % Area Covered
- ✓ Connected Sets
- ✓ Tube Length Per Set

- **Nodes:** number of connected blobs with thickness exceeding maximum width; excluded from length and area measures
- **Total Node Area:** Total square microns of node area
- **Mean Node Area:** Total square microns of node area divided by the number of nodes
- **Node % Area Covered:** Total node area divided by total image area (width times height), times 100

# Summary Data (site-by-site measurements)

- ✓ Total Tube Length
- ✓ Mean Tube Length
- ✓ Total Tube Area
- ✓ Mean Tube Area
- ✓ Tube % Area Covered
- ✓ Average Tube Thickness
- ✓ Segments
- ✓ Branch Points
- ✓ Nodes
- ✓ Total Node Area
- ✓ Mean Node Area
- ✓ Node % Area Covered
- ✓ Connected Sets
- ✓ Tube Length Per Set

- **Connected Sets:** Number of distinct objects detected in the image not connected to one another (no path of connected pixels of tubes or nodes connects the objects); measures the overall connectivity of the growth network (a completely connected network would have just one connected set of pixels)
- **Tube Length Per Set:** Total tube length in microns divided by the number of connected sets



Together through life sciences.

[www.moleculardevices.com](http://www.moleculardevices.com)