

**Molecular
Devices**

Together through life sciences.

MetaXpress® Software: *Neurite Outgrowth Module*



Together through life sciences.

©2014 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

 **Molecular
Devices**

Neurite Outgrowth

- Identifies and measures processes connected to cell bodies
- Optional use of a nuclear stain for assisting cell body identification
- Cell bodies / outgrowths can be either fluorescent or transmitted light images



Module Settings

Configure Settings for Neurite Outgrowth - neurite out...

Neurite image: 470/535 FITC

Display result image: neurite

Adaptive Background Correction™ system

Illumination

Fluorescence Transmission

Cell bodies

Approximate max width: 30 μm = 10 pixels

Intensity above local background: 200 graylevels

Minimum area: 200 μm^2 = 20 pixels

Nuclear stain (optional)

Nuclear image: 405/46E

Display result: [None]

Approximate min width: 10 μm = 3 pixels

Approximate max width: 30 μm = 10 pixels

Intensity above local background: 100 graylevels

Outgrowths

Maximum width: 20 μm = 6 pixels

Intensity above local background: 30 graylevels

Minimum cell growth to log as significant: 300 μm = 96 pixels

Configure Summary Log... Configure Data Log (Cells)...

Save Settings... Load Settings... Set to Defaults Test Run Close

- **Neurite image**
- Select the wavelength for the **neurites** (outgrowths and cell bodies)

Module Settings

Configure Settings for Neurite Outgrowth - neurite out...

Neurite image: 470/535 FITC

Display result image: neurite

Adaptive Background Correction™ system

Illumination
 Fluorescence Transmission

Cell bodies

Approximate max width: 30 $\mu\text{m} = 10$ pixels
Intensity above local background: 200 graylevels
Minimum area: 200 $\mu\text{m}^2 = 20$ pixels

Nuclear stain (optional)

Nuclear image: 405/46E

Display result: [None]

Approximate min width: 10 $\mu\text{m} = 3$ pixels
Approximate max width: 30 $\mu\text{m} = 10$ pixels
Intensity above local background: 100 graylevels

Outgrowths

Maximum width: 20 $\mu\text{m} = 6$ pixels
Intensity above local background: 30 graylevels
Minimum cell growth to log as significant: 300 $\mu\text{m} = 96$ pixels

Configure Summary Log... Configure Data Log (Cells)...

Save Settings... Load Settings... Set to Defaults Test Run Close

- Display result image
- Leave “Display result image” deselected (this is generally only used when journaling)

Module Settings

Configure Settings for Neurite Outgrowth - neurite out...

Neurite image: 470/535 FITC

Display result image: neurite

Adaptive Background Correction™ system

Illumination

Fluorescence Transmission

Cell Bodies

Approximate max width: 30 $\mu\text{m} = 10$ pixels

Intensity above local background: 200 graylevels

Minimum area: 200 $\mu\text{m}^2 = 20$ pixels

Nuclear stain (optional)

Nuclear image: 405/46E

Display result: [None]

Approximate min width: 10 $\mu\text{m} = 3$ pixels

Approximate max width: 30 $\mu\text{m} = 10$ pixels

Intensity above local background: 100 graylevels

Outgrowths

Maximum width: 20 $\mu\text{m} = 6$ pixels

Intensity above local background: 30 graylevels


Minimum cell growth to log as significant: 300 $\mu\text{m} = 96$ pixels

Configure Summary Log... Configure Data Log (Cells)...

Save Settings... Load Settings... Set to Defaults Test Run Close


- **Illumination**
- Select **Fluorescence** or **Transmission** as appropriate for the neurite image
- **Fluorescence**: light objects on dark background
- **Transmission** (Brightfield or Phase Contrast): dark objects on light background

Module Settings




- **Cell bodies**
- **Set the Approximate max width** for the cell bodies that you want to detect
- The width is the short axis of a cell body (in μm)
- The region tools can be used to measure widths
- If the setting is too large, clusters of cell bodies may be joined together
- If the setting is too small, thicker outgrowths or branch points may be identified as cell bodies


Module Settings



- **Cell bodies**
- The **intensity above local background** is used for finding the cell bodies
- This value is a minimum and should be set slightly lower than the difference in intensity between a dim cell and its local background
- Draw a line across a cell into the background and use Measure > Linescan to determine intensity values; or simply mouse over the cell and the background and view the intensity values



Module Settings



- **Cell bodies**
- The **minimum area** is an additional criteria to separate cell bodies from branch points / thick outgrowths
- Use the elliptical region tool to measure the area of a small cell in the image
- If the area is set too large, smaller cell bodies will be missed

Module Settings

Configure Settings for Neurite Outgrowth - neurite out...

Neurite image: 470/535 FITC Adaptive Background Correction™ system

Display result image: neurite

Illumination
 Fluorescence Transmission

Cell bodies

Approximate max width: 30 $\mu\text{m} = 10$ pixels
Intensity above local background: 200 graylevels
Minimum area: 200 $\mu\text{m}^2 = 20$ pixels

Nuclear stain (optional)

Nuclear image: 405/46E
 Display result: [None]

Approximate min width: 10 $\mu\text{m} = 3$ pixels
Approximate max width: 30 $\mu\text{m} = 10$ pixels
Intensity above local background: 100 graylevels


Outgrowths

Maximum width: 20 $\mu\text{m} = 6$ pixels
Intensity above local background: 30 graylevels
Minimum cell growth to log as significant: 300 $\mu\text{m} = 96$ pixels

Configure Summary Log... Configure Data Log (Cells)...

Save Settings... Load Settings... Set to Defaults Test Run Close

- **Outgrowths**
- The **maximum width** is used to help distinguish the beginning of outgrowths from cell bodies or branch points
- Use the line region tool to measure the width of a thick outgrowth



Module Settings

Configure Settings for Neurite Outgrowth - neurite out...

Neurite image: 470/535 FITC

Display result image: neurite

Adaptive Background Correction™ system

Illumination

Fluorescence Transmission

Cell bodies

Approximate max width: 30 $\mu\text{m} = 10$ pixels

Intensity above local background: 200 graylevels

Minimum area: 200 $\mu\text{m}^2 = 20$ pixels

Nuclear stain (optional)

Nuclear image: 405/46E

Display result: [None]

Approximate min width: 10 $\mu\text{m} = 3$ pixels

Approximate max width: 30 $\mu\text{m} = 10$ pixels

Intensity above local background: 100 graylevels

Outgrowths

Maximum width: 20 $\mu\text{m} = 6$ pixels

Intensity above local background: 30 graylevels


Minimum cell growth to log as significant: 300 $\mu\text{m} = 96$ pixels

Configure Summary Log... Configure Data Log (Cells)...

Save Settings... Load Settings... Set to Defaults Test Run Close

- **Outgrowths**
- The **intensity above local background** is used for finding the outgrowths
- This value is a minimum and should be set slightly lower than the difference in intensity between the dimmest part of the outgrowth (typically the end) and its local background
- Measure the intensity for outgrowths just as for the cell bodies

Module Settings



- **Outgrowths**
- The **minimum cell growth to log as significant** is used for scoring cells as positive or negative for outgrowth
- This setting does not affect the detection of the outgrowths
- This is a length measurement and can be determined using the region tools

Module Settings

Configure Settings for Neurite Outgrowth - neurite out...

Neurite image: 470/535 FITC

Display result image: neurite

Adaptive Background Correction™ system

Illumination

Fluorescence Transmission

Cell bodies

Approximate max width: 30 μm = 10 pixels

Intensity above local background: 200 graylevels

Minimum area: 200 μm² = 20 pixels

Nuclear stain (optional)

Nuclear image: 405/465 DAPI

Display result: [None]

Approximate min width: 10 μm = 3 pixels

Approximate max width: 30 μm = 10 pixels

Intensity above local background: 100 graylevels

Outgrowths

Maximum width: 20 μm = 6 pixels

Intensity above local background: 30 graylevels

Minimum cell growth to log as significant: 300 μm = 96 pixels

Configure Summary Log... Configure Data Log (Cells)...

Save Settings... Load Settings... Set to Defaults Test Run Close

- **Nuclear stain (optional)**
- If using this option, select the wavelength for the **nuclei**
- Using the nuclear stain can improve cell body identification in certain cases (dense sample, no non-neuronal cells)
- If there are non-neuronal cells present, the nuclear stain may make cell body identification worse
- For a new assay, it is recommended to test the analysis with and without the nuclear stain option

Module Settings

Configure Settings for Neurite Outgrowth - neurite out...

Neurite image: 470/535 FITC

Display result image: neurite

Adaptive Background Correction™ system

Illumination

Fluorescence Transmission

Cell bodies

Approximate max width: 30 $\mu\text{m} = 10$ pixels

Intensity above local background: 200 graylevels

Minimum area: 200 $\mu\text{m}^2 = 20$ pixels

Nuclear stain (optional)

Nuclear image: 405/465 DAPI

Display result: [None]

Approximate min width: 10 $\mu\text{m} = 3$ pixels

Approximate max width: 30 $\mu\text{m} = 10$ pixels

Intensity above local background: 100 graylevels

Outgrowths

Maximum width: 20 $\mu\text{m} = 6$ pixels

Intensity above local background: 30 graylevels


Minimum cell growth to log as significant: 300 $\mu\text{m} = 96$ pixels

Configure Summary Log... Configure Data Log (Cells)...

Save Settings... Load Settings... Set to Defaults Test Run Close

- Nuclear stain (optional)
- Leave “**Display result**” deselected (this is generally only used when journaling)


Module Settings



- **Nuclear stain (optional)**
- **Set the Approximate min width and Approximate max width** for the range of nuclei that you want to detect
- The width is the short axis of a nucleus (in μm)
- The region tools can be used to measure widths
- Much smaller cells will be ignored
- Much larger cells will be split

Module Settings


Effects of varying nuclear width settings



Min width too small: splits nuclei



Min width too large: omits smaller nuclei



Max width too small: may shrink nuclear boundaries



Max width too large: may slightly enlarge nuclear boundaries

Module Settings

Configure Settings for Neurite Outgrowth - neurite out... Adaptive Background Correction™ system

Neurite image: 470/535 FITC

Display result image: neurite

Illumination
 Fluorescence Transmission

Cell bodies

Approximate max width: 30 μm = 10 pixels

Intensity above local background: 200 graylevels

Minimum area: 200 μm² = 20 pixels

Nuclear stain (optional)

Nuclear image: 405/465 DAPI

Display result: [None]

Approximate min width: 10 μm = 3 pixels

Approximate max width: 30 μm = 10 pixels

Intensity above local background: 100 graylevels

Outgrowths

Maximum width: 20 μm = 6 pixels

Intensity above local background: 30 graylevels

Minimum cell growth to log as significant: 300 μm = 96 pixels


Configure Summary Log... Configure Data Log (Cells)...

Save Settings... Load Settings... Set to Defaults Test Run Close

- Nuclear stain (optional)
- The intensity above local background is used for finding the nuclei
- This value is a minimum and should be set slightly lower than the difference in intensity between a dim nucleus and its local background
- Measure the nuclear intensity just as for the cell bodies

Regions for Measurement

Processes (4)
connected to
cell body



Summary Data (site-by-site measurements)

- ✓ Number of Cells
- ✓ Total Outgrowth
- ✓ Mean Outgrowth Per Cell
- ✓ Total Processes
- ✓ Mean Processes Per Cell
- ✓ Total Branches
- ✓ Mean Branches Per Cell
- ✓ Total Cell Body Area
- ✓ Mean Cell Body Area
- ✓ Straightness
- ✓ Cells Significant Growth
- ✓ %Cells Significant Growth
- ✓ Mean Outgrowth Average Intensity

- **Number of Cells:** Number of cell bodies in the image
- **Total Outgrowth:** Total length of skeletonized outgrowth in mm (corrected for diagonal lengths)
- **Mean Outgrowth Per Cell:** Average skeletonized outgrowth in mm corrected for diagonal lengths divided by the number of cells.

Summary Data (site-by-site measurements)

- ✓ Number of Cells
- ✓ Total Outgrowth
- ✓ Mean Outgrowth Per Cell
- ✓ Total Processes
- ✓ Mean Processes Per Cell
- ✓ Total Branches
- ✓ Mean Branches Per Cell
- ✓ Total Cell Body Area
- ✓ Mean Cell Body Area
- ✓ Straightness
- ✓ Cells Significant Growth
- ✓ %Cells Significant Growth
- ✓ Mean Outgrowth Average Intensity

- **Total Processes:** Number of outgrowths in the image that are connected to cell bodies
- **Mean Processes Per Cell:** Total Processes divided by Number of Cells
- **Total Branches:** Total number of branching junctions in the image
- **Mean Branches Per Cell:** Total Branches divided by Number of Cells

Summary Data (site-by-site measurements)

- ✓ Number of Cells
- ✓ Total Outgrowth
- ✓ Mean Outgrowth Per Cell
- ✓ Total Processes
- ✓ Mean Processes Per Cell
- ✓ Total Branches
- ✓ Mean Branches Per Cell
- ✓ Total Cell Body Area
- ✓ Mean Cell Body Area
- ✓ Straightness
- ✓ Cells Significant Growth
- ✓ %Cells Significant Growth
- ✓ Mean Outgrowth Average Intensity

- **Total Cell Body Area:** Total μm^2 of the cell bodies in the image (excluding outgrowths)
- **Mean Cell Body Area:** Total Cell Body Area divided by the Number of Cells
- **Straightness:** Ratio varying between 0 (not straight) and 1 (perfectly straight) defined as end-to-end Euclidean distance between segment junctions divided by corresponding actual neurite curve length (the sum of end-to-end lengths divided by the sum of curve lengths), averaged over all of the cells in the image

Summary Data (site-by-site measurements)

- ✓ Number of Cells
- ✓ Total Outgrowth
- ✓ Mean Outgrowth Per Cell
- ✓ Total Processes
- ✓ Mean Processes Per Cell
- ✓ Total Branches
- ✓ Mean Branches Per Cell
- ✓ Total Cell Body Area
- ✓ Mean Cell Body Area
- ✓ Straightness

- ✓ Cells Significant Growth
- ✓ %Cells Significant Growth
- ✓ Mean Outgrowth Average Intensity

- **Cells Significant Growth:** Number of cells in the image with outgrowth greater than the threshold length specified in the settings
- **%Cells Significant Growth:** Cells Significant Growth, divided by the Number of Cells, times 100
- **Mean Outgrowth Average Intensity:** The average pixel intensity of the neurite stain over all the outgrowths detected in the image

Cell Data (cell-by-cell measurements)

- ✓ Cell: Assigned Label #
- ✓ Cell: Total Outgrowth
- ✓ Cell: Processes
- ✓ Cell: Mean Process Length
- ✓ Cell: Median Process Length
- ✓ Cell: Max Process Length
- ✓ Cell: Branches
- ✓ Cell: Straightness
- ✓ Cell: Cell Body Area
- ✓ Cell: Mean Outgrowth Intensity

- **Cell: Assigned Label #** – Cell label number (1 through total cell number), corresponds to intensity value in result image (if used)
- **Cell: Total Outgrowth** – Total amount of skeletonized outgrowth in μm (corrected for diagonal lengths) associated with the cell
- **Cell: Processes** – Number of outgrowths that connect to the cell body

Cell Data (cell-by-cell measurements)

- ✓ Cell: Assigned Label #
- ✓ Cell: Total Outgrowth
- ✓ Cell: Processes
- ✓ Cell: Mean Process Length
- ✓ Cell: Median Process Length
- ✓ Cell: Max Process Length
- ✓ Cell: Branches
- ✓ Cell: Straightness
- ✓ Cell: Cell Body Area
- ✓ Cell: Mean Outgrowth Intensity

- **Cell: Mean Process Length** – Total outgrowth (in μm) divided by number of processes of the cell
- **Cell: Median Process Length** – Median value of the outgrowth lengths (in μm) associated with the cell's various processes
- **Cell: Max Process Length** – Maximum value of the outgrowth lengths (in μm) associated with the cell's various processes

Cell Data (cell-by-cell measurements)

- ✓ Cell: Assigned Label #
- ✓ Cell: Total Outgrowth
- ✓ Cell: Processes
- ✓ Cell: Mean Process Length
- ✓ Cell: Median Process Length
- ✓ Cell: Max Process Length
- ✓ Cell: Branches
- ✓ Cell: Straightness
- ✓ Cell: Cell Body Area
- ✓ Cell: Mean Outgrowth Intensity

- **Cell: Branches** – Number of branching junctions of all the processes connected to the cell.
- **Cell: Straightness** – Ratio varying between 0 (not straight) and 1 (perfectly straight) defined as end-to-end Euclidean distance between the cell's segment junctions divided by corresponding actual neurite curve length (the sum of end-to-end lengths divided by the sum of curve lengths)

Cell Data (cell-by-cell measurements)

- ✓ Cell: Assigned Label #
- ✓ Cell: Total Outgrowth
- ✓ Cell: Processes
- ✓ Cell: Mean Process Length
- ✓ Cell: Median Process Length
- ✓ Cell: Max Process Length
- ✓ Cell: Branches
- ✓ Cell: Straightness
- ✓ Cell: Cell Body Area
- ✓ Cell: Mean Outgrowth Intensity

- **Cell: Cell Body Area** – Total area in μm^2 of the cell body (excluding outgrowths)
- **Cell: Mean Outgrowth Intensity** - The average pixel intensity of the neurite stain over all the outgrowths for this cell



Together through life sciences.

www.moleculardevices.com