



Swinging-Bucket Rotor and Buckets, Titanium

For use in instruments classified: **H R S**

Major applications: Separation of subcellular particles and viruses in density gradients.

Max. RPM	Max. <i>g</i>	<i>k</i> Factor	Number of Tubes Volume/Size	Rotor Capacity
32 000	175 000	204	6 x 38.5 mL 1 x 3.5 in 25 x 89 mm	231 mL

No. 369694. SW 32 Ti Rotor Package. Includes Rotor, 1 box 326823 Polyallomer Tubes, 2 boxes 344058 Ultra-Clear Tubes, 1 box 355631 Polycarbonate Tubes, 1 each 331186 Bucket Holder Rack, 306812 Spinkote Lubricant, 335148 Vacuum Grease, 335456 Overspeed Disk, 978354 hooked O-ring Removal Tool, and 24 each 812715 O-rings.

No. 369650. SW 32 Ti Rotor Assembly.

Tubes and Bottles

Tube Style/Material	Vol. (mL) <i>g-Max</i>	Part No.	Quantity	Size (mm)	Floating Spacers (qty. 1)	Adapters (qty. 6)	<i>g</i> -Force	<i>k</i> Factor	Max. Speed
OptiSeal Polyallomer	32.4	361625	56	25 x 77	392833 ^{§§}	—	175 000	155	32 000
Quick-Seal Polyallomer Bell-Top	15.0 27.0 33.5	343664 343665 344623	50 50 50	25 x 38 25 x 64 25 x 83	355536 355536 355536	— — —	175 000 175 000 175 000	81 146 192	32 000 32 000 32 000
Quick-Seal <i>k</i> onical Polyallomer Bell-Top	8.4 22.5 28.0	358652 358654 358651	50 50 50	25 x 38 25 x 76 25 x 83	355536 355536 355536	358156 358156 358156	175 000 175 000 175 000	74 155 198	32 000 32 000 32 000
Thinwall <i>k</i> onical Polyallomer	25.5 31.5	358125 358126	50 50	25 x 76 25 x 89	— —	358156 358156	175 000 175 000	156 156	32 000 32 000
Thickwall Polyallomer	31.0	355642	25	25 x 89	—	—	175 000	198	32 000
Thinwall Polyallomer	38.5	326823	50	25 x 89	—	—	175 000	204	32 000
Thickwall Polycarbonate	31.0	355631	25	25 x 89	—	—	175 000	202	32 000
Ultra-Clear	38.5	344058	50	25 x 89	—	—	175 000	204	32 000

Note: To help you locate *g-Max* tubes in the rotor listings, these tubes are highlighted with bold, italic type.

^{§§} Set of 2.

Rotor Replacement Parts

331186 Bucket Holder Rack
332400 Rotor Stand
335456 Overspeed Disk, 32 000 rpm
369643 Bucket Cap (each)
369692 Bucket Assembly, with Caps and O-rings, set of 6
812715 O-ring for Bucket

Extra Buckets for the SW 32 Ti

The buckets of the SW 32.1 can be used with the SW 32 rotor head.
369693 Bucket Assembly, with Caps and O-rings, set of 6

Adapters/Spacers

355536

358156

392833

