


1



**Swinging-Bucket Rotor, Titanium Head, and Buckets (black buckets)**

For use in instruments classified: **H R S**

Major applications: Rate-zonal and isopycnic centrifugation of viruses, rate-zonal centrifugation of RNA.

Max. RPM	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
41 000	288 000	124	6 x 13.2 mL 1/16 x 3 1/2 in 14 x 89 mm	79.2 mL

**No. 331336.** SW 41 Ti Rotor Package. Includes Rotor, 2 boxes 331372 Polyallomer Tubes, 2 boxes 344059 Ultra-Clear Tubes, 1 each 330070 Hinge Pin Tool, 331313 Bucket Holder Rack, 306812 Spinkote Lubricant, 335148 Vacuum Grease, 330335 Overspeed Disk, and 24 each 331309 Gaskets for Buckets.

**No. 331362.** SW 41 Ti Rotor Assembly.

**Tubes and Bottles**

Tube Style/Material	Vol. (mL) <i>g-Max</i>	Part No.	Quantity	Size (mm)	Floating Spacers (qty. 1)	Adapters (qty. 6)	g-Force	k Factor	Max. Speed
Quick-Seal	<b>3.5</b>	355870	50	14 x 25	<b>355534</b>	—	288 000	27	41 000
Polyallomer	<b>5.9</b>	355537	50	14 x 47	<b>355534</b>	—	288 000	55	41 000
Quick-Seal <i>k</i> onical	<b>4.0</b>	358650	50	14 x 48	<b>355534</b>	<b>358154*</b>	288 000	57	41 000
Polyallomer	8.0	358649	50	14 x 89	355534	358154*	288 000	124	41 000
Thinwall Polyallomer	13.2	331372	50	14 x 89	—	—	288 000	124	41 000
Thinwall <i>k</i> onical Polyallomer	10.0	358120	50	14 x 89	—	358154*	284 000	124	41 000
Ultra-Clear	13.2	344059	50	14 x 89	—	—	288 000	124	41 000

*Note: To help you locate g-Max tubes in the rotor listings, these tubes are highlighted with bold, italic type.*

\* Recommend P/N 354468 Extractor Tool (konical tube adapters).

**Rotor Replacement Parts**

- 330335 Overspeed Disk, 41 000 rpm
- 331309 Gasket for Bucket
- 331763 Caps, matched set of 6
- 331313 Bucket Holder Rack
- 332400 Rotor Stand
- 333790 Bucket Assembly, Black Titanium, with Caps and Gaskets, matched set of 6

**Adapters/Spacers**

355534

358154

